


cellenCHIP 384 - 3'RNA-seq Kit Sequencing Guidelines



Copyright notice

© Copyright 2022 by Cellenion SASU. All rights reserved.

Information in this document is subject to change without notice and does not represent a commitment on the part of Cellenion SASU. Every effort has been made to supply complete and accurate information; however, Cellenion SASU assumes no responsibility and will not be liable for any error, omissions, damage, or loss that might result from any use of this manual or the information contained therein. Any information given may not be used in any other form. No part of this manual may be copied or distributed, transmitted, transcribed, stored in a retrieval system or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, manual or otherwise, or disclosed to third parties without the expressed written permission from:

CELLENION SASU
BioSerra 2
60 Avenue Rockefeller
69008 Lyon
France
ticket@cellenion.com

Disclaimers

The information contained in this document is the proprietary and exclusive property of Cellenion SASU except as otherwise indicated. No part of this document, in whole or in part, may be reproduced, stored, transmitted, or used for design purposes without the prior written permission of Cellenion SASU.

The information in this document is provided for informational purposes only. Cellenion SASU specifically disclaims all warranties, express or limited, including, but not limited, to the implied warranties of merchantability and fitness for a particular purpose, except as provided for in a separate software license agreement.

Trademarks

Cellenion®, cellenCHIP™

Version History

| REVISION CHART | | | |
|----------------|--|-----------------------------------|----------------|
| Version | Author(s) | Description of Version | Date Completed |
| 0.5 | J.W. Bagnoli | First creation of the document | 23-May-2022 |
| 0.9 | J.W. Bagnoli F. Monjaret | Updated figures, tables, and text | 16.06.2022 |
| 1.0 | J.W. Bagnoli F. Monjaret S. Ruiz L. Jarassier | First document release | 04.07.2022 |
| 1.1 | J.W. Bagnoli | Updated figures | 28.07.2022 |

Table of Contents

| | |
|---|----|
| I. Library Format | 3 |
| 1. Library format generated with the “cellenCHIP 384 – 3’RNA-seq Kit”... | 3 |
| 2. Compatible sequencing platforms..... | 3 |
| 3. Reads and indices of the library..... | 4 |
| II. Sequencing set up | 4 |
| 1. Single indexed sequencing scenario (recommended) | 4 |
| a. “cellenCHIP 384 – 3’RNA-seq Kit” generated libraries carry only the i7 index..... | 4 |
| b. Detailed sequencing workflow..... | 4 |
| c. Special attention is required for Read 1: 18 or 26 cycles? | 6 |
| 2. Dual indexed sequencing scenario (compatible but not optimal) | 6 |
| a. Dual-indexed sequencing is compatible with “cellenCHIP 384 – 3’RNA-seq Kit” generated libraries..... | 6 |
| b. Detailed dual indexed sequencing workflow using both workflow A or B from Illumina..... | 6 |
| 3. Recommended sequencing set up..... | 9 |
| III. Index and barcode sequences | 9 |
| 1. Index sequences..... | 9 |
| 2. Barcodes sequences..... | 10 |
| IV. Sequencing depth | 15 |

I. Library format

1. Library format generated with the “cellenCHIP 384 – 3'RNA-seq Kit”

A detailed overview of the fragments generated by the “cellenCHIP 384 – 3'RNA-seq Kit” is presented in [Figure 1](#) as well as [Table 1](#).



Figure 1: Graphical scheme of “cellenCHIP 384 – 3'RNA-seq Kit” derived libraries.

Table 1: Description of sequence blocks within final library derived from the “cellenCHIP 384 – 3'RNA-seq Kit”

| Name | Sequence (5'-3') | Length | Description |
|-------|------------------------------------|---------------|--|
| P7 | CAAGCAGAAGACGGCATACGAGAT | 24 bp | Illumina P7 adapter for flowcell binding |
| i7 | Variable, known | 8 bp | i7 index, discriminates different cellenCHIPs 384 |
| Read2 | GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG | 34 bp | Nextera Read 2 read start sequence |
| cDNA | Variable | -50 to~900 bp | cDNA fragment derived from mRNA molecule |
| PolyA | AAAAAAAAAAAAAAAAAAAAAAAAAAAA | ~30 bp | 3' polyadenylation sequence of mRNA/oligodT sequence of oligodT primers, can slightly vary in length |
| UMI | Variable, unknown | 8 bp | Unique molecular identifier |
| BC | Variable, known | 10 bp | Cell barcode, discriminating each cell on a cellenCHIP 384 |
| Read1 | AGATCGGAAGAGCGTCGTGTAGGAAAGA | 29 bp | TruSeq Read 1 read start sequence |
| P5 | GTGTAGATCTCGGTGGTCGCCGTATCATT | 29 bp | Illumina P5 adapter for flowcell binding |

The library is designed for single indexed sequencing runs, but is also compatible with dual indexed sequencing runs (see [section II.2.](#)).

2. Compatible sequencing platforms

The “cellenCHIP 384 – 3'RNA-seq Kit” generates sequencing ready libraries uniquely suited for Illumina’s next generation “Sequencing by Synthesis” (SBS) platforms. They are compatible with the Illumina iSeq, MiniSeq, MiSeq, NextSeq, HiSeq and NovaSeq systems and do not require custom sequencing primers. Other sequencing methods (e.g. BGI/MGI, IonTorrent, 454 Pyrosequencing, SOLID sequencing, Nanopore sequencing...) are not supported. To sequence on such platforms please refer to the corresponding potential conversion kits provided by the relevant suppliers.

If you are not familiar with Illumina sequencing we highly recommend to have a look at Illumina’s [“Introduction to SBS Technology”](#) as well as their [“Index Sequencing Overview Guide”](#) and additionally to the detailed explanation of the group of Sarah Teichmann on [“Github”](#).

3. Reads and indices of the library

- **Read 1 (Illumina TruSeq Read 1)** is used to read out the cell barcode (BC), which carries the information of the cell/nanowell the fragment originated from, as well as the UMI (Unique molecular identifier), which carries the information from which original RNA molecule the fragment originated from. Please note that while the cell barcode (10 nt) and the UMI (8 nt) only require a 18 nt to be sequenced, Illumina sequencers often require a minimum of 26 cycles in Read 1. Please refer to [section II.1.c.](#) for further information.
- **Index Read 1 (i7 index)** is used to read out the cellenTAG index, which distinguishes different cellenCHIPs when present on the same flowcell and requires a standard 8 cycle index read.
- Finally, **Read 2 (Illumina Nextera Read 2)** is used to read out the sequence information of the cDNA. This read can be used for mapping towards the corresponding genome and acquire the expression data. This read should be of a minimum length of 50 cycles for well annotated genomes and should not exceed 100 cycles.

Please note that the libraries do not contain an Index 2 (i5). For further information please refer to [section II.2.](#)

II. Sequencing set up

1. Single indexed sequencing scenario (recommended)


a. “cellenCHIP 384 – 3'RNA-seq Kit” generated libraries carry only the i7 index

In a single indexed sequencing run, only index 1 (i7) is read. Libraries generated with the “cellenCHIP 384 – 3'RNA-seq Kit” carry only a single Illumina index (i7, Index 1) which discriminates between different cellenCHIPs when pooled on the same flowcell or lane. The 8 nt long i7 index, introduced during tagmentation using the cellenTAG enzyme, is read in the second read of the sequencing workflow using the Index 1 primer. For a detailed list of the different Indices integrated by the cellenTAG enzymes please refer to [Table 3](#).

b. Detailed sequencing workflow

[Figure 2](#) depicts a detailed scheme of the sequencing in a single indexed run. In short, after cluster generation, the TruSeq Read 1 sequencing primer anneals to the template and initiates the SBS reaction. Afterwards, the Read 1 product is removed and the Index 1 (i7) sequencing primer is annealed to the same template strand, producing the Index 1 (i7) Read. As before, the generated product is removed, and the original template strand is used to regenerate the complementary strand. Then, the original template strand is removed to allow the Nextera Read 2 sequencing primer to bind, subsequently followed by the SBS reaction to obtain Read 2.

A)



B)




Figure 2: Schematic representation of single indexed sequencing of a “cellenCHIP 384 – 3'RNA-seq Kit” library.

A) Schematic representation of the final library structure and the corresponding sequencing primers.

B) Schematic representation of each SBS reaction on a Illumina flowcell. Each Illumina sequencing read (Read 1, Index 1, Read 2) is depicted in a separate box.

- Barcode (BC) sequence will depend on cell/nanowell of origin. Shown in this figure: barcode sequence of barcode 1.

- UMI sequence will depend on initial RNA molecule. Shown in this figure: random sequence.

- i7 sequence will depend on cellenCHIP of origin. Shown in this figure: index sequence of cellenTAG_1.

- cDNA sequence will depend on original RNA fragment captured. Shown in this figure: 50 nt sequence of human FLT3 gene (exon 24).

c. Special attention is required for Read 1: 18 or 26 cycles?

As mentioned before, Illumina sequencers often require a minimum of 26 cycles in Read 1, as it is used by the system to calculate quality metrics correctly ([Why sequencing 26 or more cycles in Read 1 is recommended](#)). For libraries generated with the “cellenCHIP 384 – 3'RNA-seq Kit” this means that additionally to the cell barcode and UMI (18 nt in total) a further 7 nt are read. Due to the structure of the library, these nucleotides will almost exclusively be “T” (derived from the initial PolyA tail of the mRNA). While this does not affect the quality of the UMI and barcode sequencing, it can, depending on the instrument used, be problematic for the quality metrics, as highly homopolymeric sequences can cause problems with cluster identification as well as phasing correction. This is particularly a problem with older Illumina systems working with non-patterned flowcells (e.g. MiSeq, HiSeq). A common way to reduce this problem is to add highly variable sequences. This can be achieved by either pooling the libraries with other libraries that do not have homopolymeric stretches in this region or specifically by adding PhiX library to the sequencing run. The exact amount of how much of highly variable sequences is required for a good quality sequencing run depends on the type of machine. Hence, before sequencing, make sure that the operator of the sequencing machine is aware of the homopolymer stretch in Read 1 and is trained to adjust the parameters accordingly.

2. Dual indexed sequencing scenario (compatible but not optimal)

a. Dual-indexed sequencing is compatible with “cellenCHIP 384 – 3'RNA-seq Kit” generated libraries

Libraries generated with the cellenCHIP 384 – 3'RNA-seq Kit are not designed to be sequenced using a dual indexing set up, but there is no restriction to use such a set up. Detailed outcomes are presented in the paragraph below.

b. Detailed dual indexed sequencing workflow using both workflow A or B from Illumina.

Dual-indexed sequencing on a paired-end flow cell follows one of two possible workflows, depending on the system and software ([Figure 3](#) and [Figure 4](#)).

- The forward strand workflow (A) is performed on the NovaSeq 6000 with v1.0 reagent kits, MiniSeq with Rapid Reagent kits, MiSeq, HiSeq 2500, and HiSeq 2000.
- The reverse complement workflow (B) is performed on the iSeq 100, MiniSeq with Standard reagent kits, NextSeq Systems, NovaSeq 6000 with v1.5 reagent kits, HiSeq X, HiSeq 4000, and HiSeq 3000.

If required, libraries generated with the “cellenCHIP 384 – 3'RNA-seq Kit” can be sequenced on a dual indexed sequencing run. In this case, instead of reading an actual Index 2 (i5) sequence, either the TruSeq Read 1 Read start sequence (Workflow A) or the P5 adapter (Workflow B) is read. These sequences are equal in all libraries and can therefore not be used to distinguish between cellenCHIP libraries. However, as these are normally not used as Index 2 sequences, they can be used to distinguish libraries generated with the “cellenCHIP 384 – 3'RNA-seq Kit” from libraries of other sources which carry an Index 2 index sequence. Please refer to [Table 4](#) to obtain these sequences.

It should be noted that there are 2 differences between the workflows. Firstly, workflow A reads the Index 2 on the forward strand and workflow B on the reverse complement strand. Secondly while workflow A uses the grafted P5 oligo as the SBS initiating primer, workflow B uses a primer binding to the Read 1 read start sequence (see [Figure 3](#) and [Figure 4](#)).





Figure 3: Schematic representation of dual indexed sequencing of a “cellenCHIP 384 – 3’RNA-seq Kit” library using the **forward (A) workflow**.

A) Schematic representation of the final library structure and the corresponding sequencing primers.

B) Schematic representation of each SBS reaction on a Illumina Flowcell. Each Illumina sequencing read (Read 1, Index 1, Index 2, Read 2) is depicted in a separate box.

- BC sequence depend on nanowell of origin. Shown in this figure: barcode sequence of barcode 1
- UMI sequence depend on initial RNA molecule. Shown in this figure: random sequence.
- i7 sequence depend on cellenCHIP of origin. Shown in this figure: index sequence of cellenTAG _1
- i5 sequence is not a real i5 index but a part of the Read 1 read start sequence, as depicted. It is the same for every fragment.
- cDNA sequence depend on original RNA fragment captured. Shown in this figure: 50 nt sequence of human FLT3 gene (Exon 24).

A)



B)




Figure 4: Schematic representation of Dual indexed Sequencing of a “cellenCHIP 384 – 3'RNA-seq Kit” library using the **reverse complement (B) workflow**.

A) Schematic representation of the final library structure and the corresponding sequencing primers.

B) Schematic representation of each SBS reaction on a Illumina Flowcell. Each Illumina sequencing read (Read 1, Index 1, Index 2, Read 2) is depicted in a separate box.

- BC sequence depend on nanowell of origin. Shown in this figure: barcode sequence of barcode 1
- UMI sequence depend on initial RNA molecule. Shown in this figure: random sequence.
- i7 sequence depend on cellenCHIP of origin. Shown in this figure: index sequence of cellenTAG _1
- i5 sequence is not a real i5 index but a part of the P5 adapter sequence, as depicted. It is the same for every fragment.
- cDNA sequence depend on original RNA fragment captured. Shown in this figure: 50 nt sequence of human FLT3 gene (Exon 24).

3. Recommended sequencing set up

We recommend the following sequencing layouts when sequencing libraries generated by the “cellenCHIP 384 – 3'RNA-seq Kit”:

Table 2: Recommended sequencing design for “cellenCHIP 384 – 3'RNA-seq Kit” derived libraries

| | Read 1 (BC/UMI) | Index 1 (i7) | Index 2 (i5) | Read 2 (cDNA) |
|-------------|------------------|-----------------|-----------------|------------------|
| Minimum | 18 cycles | 8 cycles | 0 cycles | 50 cycles |
| Maximum | 26 cycles | 8 cycles | 0 cycles | 100 cycles |
| Recommended | 26 cycles | 8 cycles | 0 cycles | 75 cycles |

III. Index and barcode sequences

1. Index sequences

In the “cellenCHIP 384 – 3'RNA-seq Kit” i7 indices are integrated to the fragments during tagmentation using the cellenTAG enzymes. Each of these enzymes carry an unique 8 nt index sequence (1 to 16) (Table 3). The i7 index can be used to distinguish between different cellenCHIP libraries on the same flowcell/lane. It also can be used to distinguish your library from any other on the sequencing run. While most sequencing runs perform 8 cycles for the index read, sometimes this is increased to 10 cycles when using a 10 nt index (e.g. unique dual indices by IDT). In this case, a further 2 nucleotides are read after the actual cellenTAG index sequences, which are part of the P7 adaptor sequence (“AT”) (see Figure 2 A, Figure 3 A and Figure 4 A). When performing a sequencing run outside your lab make sure to provide the sequencing operator with the correct index sequence(s) as depicted below. The sequences below depict the sequence necessary to be provided in the sample sheet for accurate demultiplexing.

Table 3: Index 1 (i7) sequences provided with the “cellenCHIP 384 – 3'RNA-seq Kit”

| cellenTAG | Index 1 if 8 nt are read | Index 1 if 10 nt are read |
|-----------|--------------------------|---------------------------|
| 1 | TTGGAATG | TTGGAATGAT |
| 2 | TTAATGCG | TTAATGCGAT |
| 3 | AGCTACGT | AGCTACGTAT |
| 4 | GCCTCCTG | GCCTCCTGAT |
| 5 | GGGACAAC | GGGACAACAT |
| 6 | TATCCCAC | TATCCCACAT |
| 7 | CAACTGTG | CAACTGTGAT |
| 8 | ATGACTAG | ATGACTAGAT |
| 9 | CCATATCC | CCATATCCAT |
| 10 | GTAGTCAC | GTAGTCACAT |
| 11 | CGGAGATA | CGGAGATAAT |
| 12 | CACTCTCA | CACTCTCAAT |
| 13 | GAGTTCTC | GAGTTCTCAT |
| 14 | TCGCCAGC | TCGCCAGCAT |
| 15 | TGTGACTA | TGTGACTAAT |
| 16 | TTGACGTC | TTGACGTCAT |

As mentioned previously, when performing a dual indexed run the constant sequences of the P5 or the TruSeq Read 1 sequence can be used as “index” sequences. Table 4 lists these sequences for both Workflows with a 8 or 10 cycle index read. As before, the sequences below depict the sequence necessary to be provided in the sample sheet for accurate demultiplexing.

Table 4: Index 2 (i5) sequences required to perform demultiplexing of a dual indexed Sequencing run of “cellenCHIP 384 – 3'RNA-seq Kit” derived libraries.

| Workflow | Index 2 if 8 nt are read | Index 2 if 10 nt are read |
|--|--------------------------|---------------------------|
| Workflow A (forward) | TCTTTCCC | TCTTTCCCTA |
| Workflow B (reverse complement) | AGATCTCG | AGATCTCGGT |

2. Barcodes sequences

In the “cellenCHIP 384 – 3'RNA-seq Kit” each well of the cellenCHIP 384 RTready chip harbors a unique barcoded oligodT primer. During the reverse transcription reaction, these primers are used to target polyadenylated RNA and thereby tag each molecule within a given well with the same barcode sequence. [Figures 6-9](#) depict each array of the cellenCHIP384 as described in [Figure 5](#), with their corresponding oligodT primer (e.g. oligodT_1) and the corresponding 10 nt barcode sequence.




Figure 5: Schematic overview of arrays and the distribution of barcoded oligodT primers in of the cellenCHIP 384.

Array 1

Column

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |
|-----------|---------------------------|---------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-----------|
| 1 | oligodT_1 TAACTTGGTC | oligodT_13 TCGAACACGA | oligodT_25 TTACTTCGTG | oligodT_37 ACGAGACGTC | oligodT_49 CATAAGCGAC | oligodT_61 AACTCTCTAC | oligodT_73 GTATGCCACA | oligodT_85 GCCAACAGAC | 1 |
| 2 | oligodT_2 ATAAGCGGAG | oligodT_14 TCCAAC TGAA | oligodT_26 ATTCGGCTTA | oligodT_38 TATGGCAAGC | oligodT_50 CCTAAGCTGA | oligodT_62 GCCACCTAGT | oligodT_74 AGTATCCGAT | oligodT_86 CCTCGTTGGA | 2 |
| 3 | oligodT_3 CGACCTTAAT | oligodT_15 GAACCGCATT | oligodT_27 CGCGGTATCT | oligodT_39 ACTCCGAAGC | oligodT_51 CGATGTCTTG | oligodT_63 CTCCGCCATA | oligodT_75 AACGAAGATC | oligodT_87 AATAGGTAGG | 3 |
| 4 | oligodT_4 ATGTTGCCAC | oligodT_16 ACAACAGGCT | oligodT_28 TATCCGCGTC | oligodT_40 AATCTGCAGT | oligodT_52 GAGATGAGAG | oligodT_64 TGACACGATT | oligodT_76 CCTCCTGTGA | oligodT_88 TATGCGTAGC | 4 |
| 5 | oligodT_5 CAGCAGGTCA | oligodT_17 TTACGCCAT | oligodT_29 AGGATGTCCA | oligodT_41 AGTTCAGCGC | oligodT_53 TCACCATAAG | oligodT_65 GTGCTCTAGG | oligodT_77 TAGGCCAATG | oligodT_89 CCAGTCGTCA | 5 |
| 6 | oligodT_6 TCTAGCTGCG | oligodT_18 TTCTCTTAGC | oligodT_30 GCCAGCTATC | oligodT_42 CGCCAAGGTA | oligodT_54 ACTTCCGGCA | oligodT_66 GCAATCTGCT | oligodT_78 AGTAAGGAGA | oligodT_90 AGAAGACCTA | 6 |
| 7 | oligodT_7 AGATTAGCGT | oligodT_19 TCGCGACACT | oligodT_31 CTTACGAACT | oligodT_43 CCTGACATCT | oligodT_55 AAGCGCTTAG | oligodT_67 CTCTATTGAG | oligodT_79 CCGAACAGAG | oligodT_91 CCGTGAGTTC | 7 |
| 8 | oligodT_8 AGGTACAAC | oligodT_20 TGGAGATCGA | oligodT_32 CAAGACTCAA | oligodT_44 CGAGCGTAGT | oligodT_56 ACAACACAG | oligodT_68 CACTTACGGC | oligodT_80 AGAACTAACG | oligodT_92 TCAATTCGGA | 8 |
| 9 | oligodT_9 CGCACTTCGT | oligodT_21 CCAGACCGAT | oligodT_33 TCGACGGCAA | oligodT_45 CTGGATTACT | oligodT_57 AGACTAGATG | oligodT_69 TCAGCATCTA | oligodT_81 GAGCAGCTTG | oligodT_93 CGGAACACAA | 9 |
| 10 | oligodT_10 CTTGCTATTTC | oligodT_22 CCACAGCCTT | oligodT_34 TCAGTTCGAT | oligodT_46 AGCAATTGAC | oligodT_58 CCTCTATAGA | oligodT_70 GTTGCATCAA | oligodT_82 CCGATAGTAG | oligodT_94 CTATAACGGA | 10 |
| 11 | oligodT_11 TCTCACCAGG | oligodT_23 ACCGAATAGA | oligodT_35 CTAACCTCAA | oligodT_47 ATATTCTCGG | oligodT_59 TACGCTACTT | oligodT_71 ATGACAGCAC | oligodT_83 GTGTCTACAG | oligodT_95 AGATGTAAGG | 11 |
| 12 | oligodT_12 TTCCTCACGT | oligodT_24 CGGACGATTTC | oligodT_36 GACCTACTAT | oligodT_48 GACTCATCCA | oligodT_60 CGCACTGGTA | oligodT_72 TAGCGACGGA | oligodT_84 GTACAGACCG | oligodT_96 ATGCTTACTG | 12 |
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |

Figure 6: Schematic representation of the distribution of barcoded oligodT primers in Array 1 of the cellenCHIP 384.

Array 2

Column

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |
|-----------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|----------------------------|---------------------------|-----------|
| 1 | oligodT_97 GC CGCGTTAA | oligodT_109 TAGGAAGCGG | oligodT_121 ACGGTCGCAT | oligodT_133 CCTATGTGTA | oligodT_145 TTGGATTGTC | oligodT_157 GCCGTGAACA | oligodT_169 TAAGATCGGA | oligodT_181 ATCCTAGAAC | 1 |
| 2 | oligodT_98 TTCCTTGAGG | oligodT_110 GCATCCTACT | oligodT_122 TTACACACGG | oligodT_134 CAACGTTGAC | oligodT_146 AATTGGCACA | oligodT_158 AGGAGACCAC | oligodT_170 TTAGTTGCTG | oligodT_182 GCCTAGAGTC | 2 |
| 3 | oligodT_99 GAGTCGCTTC | oligodT_111 AGTCCAGCGT | oligodT_123 CGAGATTAGA | oligodT_135 TGACATAGGT | oligodT_147 AGGCAATCGC | oligodT_159 GAAGCAAGCG | oligodT_171 AGTCAGATGC | oligodT_183 CTTACTCTCT | 3 |
| 4 | oligodT_100 GTACGCCTAA | oligodT_112 AGGCAGCCAT | oligodT_124 CAAGGATCGA | oligodT_136 TTACACCAAC | oligodT_148 TGCAACACTA | oligodT_160 AACGCGACTG | oligodT_172 AGTACAACCTG | oligodT_184 CCGACAACCT | 4 |
| 5 | oligodT_101 TAGAAGATCG | oligodT_113 GCGACTCGTA | oligodT_125 ACAACACACA | oligodT_137 CCTGGAATTA | oligodT_149 GAATCGCCTG | oligodT_161 TTCGCGGACT | oligodT_173 CTTAGTACGC | oligodT_185 GATCTGTATG | 5 |
| 6 | oligodT_102 AACCTAGTGC | oligodT_114 GTATCTGAGG | oligodT_126 AAGAGCTCGC | oligodT_138 GCAACATCAA | oligodT_150 TCGCATGACA | oligodT_162 TCGTGACGAC | oligodT_174 GTACGCATCG | oligodT_186 ATAGGCAGTG | 6 |
| 7 | oligodT_103 CGTGTATGTC | oligodT_115 CATCACCTAA | oligodT_127 ACAGAGGCGT | oligodT_139 ACAGCCAGGT | oligodT_151 TCAATACGTG | oligodT_163 TCAATCTAGC | oligodT_175 ACTTCGTTTC | oligodT_187 GAATAGAGTG | 7 |
| 8 | oligodT_104 TTCAGATCCA | oligodT_116 CGGCATATAC | oligodT_128 AACTGTACGG | oligodT_140 TAGTCCATCT | oligodT_152 ATTCTCCAGC | oligodT_164 GCGGAGTCAA | oligodT_176 CCACCTGCCT | oligodT_188 GTAACAGGAA | 8 |
| 9 | oligodT_105 CTCACCAGTT | oligodT_117 CACGAGTGTC | oligodT_129 GCGAAGTAGG | oligodT_141 GCCTCTCAGA | oligodT_153 TGTTACAGAG | oligodT_165 GATGCTAGAA | oligodT_177 TGGTGGTTGT | oligodT_189 ATGGTTGTTG | 9 |
| 10 | oligodT_106 ACTAGTAGTC | oligodT_118 CGGAGACTAT | oligodT_130 GTGAAGACAC | oligodT_142 TAGGTTCAAG | oligodT_154 TTACGGACGC | oligodT_166 ACGATGGTTA | oligodT_178 ACGCAACACA | oligodT_190 ACAATAGCAC | 10 |
| 11 | oligodT_107 AATAGACTGC | oligodT_119 AAGGTGTTAG | oligodT_131 GACTTAGGTC | oligodT_143 AAGACTCTCG | oligodT_155 CCAATGTTGA | oligodT_167 AGAGCAGTAT | oligodT_179 CCAGTGCACT | oligodT_191 CTACAAGTTC | 11 |
| 12 | oligodT_108 ATGATCAACG | oligodT_120 GCGATGGAAG | oligodT_132 AGGACGGACA | oligodT_144 GAAGATCGAC | oligodT_156 CGTACTTAGA | oligodT_168 TTGTGATCCT | oligodT_180 CATCAGGCGT | oligodT_192 AACCTGATCT | 12 |
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |

Figure 7: Schematic representation of the distribution of barcoded oligodT primers in Array 2 of the cellenCHIP 384.

Array 3

Column

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |
|-----------|----------------------------|---------------------------|----------------------------|---------------------------|----------------------------|----------------------------|---------------------------|---------------------------|-----------|
| 1 | oligodT_192 AACCTGATCT | oligodT_205 GCGCCATATA | oligodT_217 GCATAAGATC | oligodT_229 GTTCTCATCA | oligodT_241 CTTGTTACTG | oligodT_253 GATAGATCAG | oligodT_265 CGCACATTAC | oligodT_277 CAGACTGCAG | 1 |
| 2 | oligodT_193 TAGCGAGGCT | oligodT_206 CACTCGCATA | oligodT_218 AAGGTGTA CT | oligodT_230 CAATGCGTGT | oligodT_242 TCACCGTGGT | oligodT_254 GCAGAAGTCT | oligodT_266 ATCTAACCTG | oligodT_278 GTTAGTGCCT | 2 |
| 3 | oligodT_194 CCAAC TGTAT | oligodT_207 AAGGTCAATG | oligodT_219 CTTCTAACCG | oligodT_231 TGAGTGTACA | oligodT_243 CAACTACCGC | oligodT_255 TATACGATGC | oligodT_267 GCGACAATAT | oligodT_279 CGCAATAAGA | 3 |
| 4 | oligodT_195 GAGTAGTGAC | oligodT_208 GTGACGCACA | oligodT_220 TCGACCAGCA | oligodT_232 TGCTTCCAAG | oligodT_244 ACTCGTCCAG | oligodT_256 AACGCTCTGA | oligodT_268 TTGCCGTATC | oligodT_280 TCGGTCTGTT | 4 |
| 5 | oligodT_196 CCAGTATTCT | oligodT_209 TGTCAGGCTT | oligodT_221 GTGACACGTC | oligodT_233 CCTACTCTGC | oligodT_245 TCCGGAGAAC | oligodT_257 GATTGGCGAT | oligodT_269 GATCCACTTA | oligodT_281 CTAGACGAGC | 5 |
| 6 | oligodT_197 ATCATA CCGC | oligodT_210 CACGTGATGG | oligodT_222 AATGGACGCG | oligodT_234 TCGCTGAGAG | oligodT_246 TTGTGAGAAG | oligodT_258 CAAGACGGAC | oligodT_270 TCCGAGATGC | oligodT_282 GTAATACTCG | 6 |
| 7 | oligodT_198 ATTCGCACTA | oligodT_211 CCATCACGTG | oligodT_223 TCTCCGTACT | oligodT_235 CGCTAATTGC | oligodT_247 CCGAAGGTTT | oligodT_259 CCACTCAAGG | oligodT_271 ATAGAGTTCT | oligodT_283 GATTCTTAGC | 7 |
| 8 | oligodT_199 CTCTTCAGCA | oligodT_212 ACCAGGAATT | oligodT_224 TGGAAGCACA | oligodT_236 TCAGCTGTCA | oligodT_248 ATACACACCT | oligodT_260 AGCCACTCGT | oligodT_272 GCATTCCTTC | oligodT_284 TGTGATCGCG | 8 |
| 9 | oligodT_200 CTGTATGCCA | oligodT_213 TAAGGCTCAC | oligodT_225 GATAGACAAC | oligodT_237 AGACGACAGA | oligodT_249 AGCCTCTGCT | oligodT_261 GCGAAGAATT | oligodT_273 AATCGTGCGG | oligodT_285 ATCGAGAGAC | 9 |
| 10 | oligodT_201 TCTGACTATG | oligodT_214 TGGCACTGAT | oligodT_226 AGCGTCTCAT | oligodT_238 TAGGTCCAGA | oligodT_250 CGGATCATT A | oligodT_262 TGTGCGGTGCT | oligodT_274 TACCGGAGCA | oligodT_286 CTGAACGTAT | 10 |
| 11 | oligodT_202 TGTGACGTAT | oligodT_215 CATATGGTGC | oligodT_227 CTCTATACTG | oligodT_239 TAACCTGAAC | oligodT_251 CTCTCTCTAG | oligodT_263 TCTACAGCTA | oligodT_275 AGCACCAGTC | oligodT_287 CACTAGACTT | 11 |
| 12 | oligodT_203 AGAAGATTGG | oligodT_216 GCACTACCAG | oligodT_228 ATAGCCGGTT | oligodT_240 GATAGGTGAA | oligodT_252 AAGGCTCCAT | oligodT_264 CAACCATGAT | oligodT_276 TACCTAGACT | oligodT_288 CGTCGCCTAT | 12 |
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |

Figure 8: Schematic representation of the distribution of barcoded oligodT primers in Array 3 of the cellenCHIP 384.

Array 4

Column

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |
|-----------|----------------------------|----------------------------|---------------------------|---------------------------|---------------------------|----------------------------|---------------------------|---------------------------|-----------|
| 1 | oligodT_289 CAACCATACA | oligodT_301 TCGACGTTGG | oligodT_313 AATCGAACTC | oligodT_325 ACTTCAACAG | oligodT_337 CGAGTTCGCA | oligodT_349 GACAATGCCG | oligodT_361 GTGCTATCGA | oligodT_373 GTCTTCTCGA | 1 |
| 2 | oligodT_290 GTAGGCCGTT | oligodT_302 AAGAGTACGA | oligodT_314 AGGCTGGAGA | oligodT_326 TGGACAATAG | oligodT_338 TCCTACGTGC | oligodT_350 AAGTGGCATT | oligodT_362 AGAACCTACA | oligodT_374 ATACACAGAG | 2 |
| 3 | oligodT_291 CGGATTGATC | oligodT_303 CAGGAGAAGT | oligodT_315 TAAGCCACGC | oligodT_327 TCACACGGAT | oligodT_339 CCACGATAAC | oligodT_351 GCAACCTGAG | oligodT_363 ATGGTGTATC | oligodT_375 AGGTCTATGG | 3 |
| 4 | oligodT_292 ACTGGCAAGA | oligodT_304 GACAATTAGC | oligodT_316 TCAATCAGAG | oligodT_328 AGGATACGGT | oligodT_340 CTTGCGAAGA | oligodT_352 AGAACCCGGAC | oligodT_364 CTCAGAATTC | oligodT_376 ACCTCATTGA | 4 |
| 5 | oligodT_293 ACAAGAACCT | oligodT_305 TTCTACCGGC | oligodT_317 AACATAAGGC | oligodT_329 GCTACAGGTC | oligodT_341 TGGCACGAAC | oligodT_353 CACACTGTTG | oligodT_365 TATGATGCAG | oligodT_377 CGTGGCAGAA | 5 |
| 6 | oligodT_294 AGAGTATGTG | oligodT_306 GTGAGAGTAT | oligodT_318 TTGCAATGCG | oligodT_330 GCCAAGTAAC | oligodT_342 AACGACATCT | oligodT_354 CTCATTCTGC | oligodT_366 GTGGAGGTAG | oligodT_378 CAACAGGTAG | 6 |
| 7 | oligodT_295 ACTTAGATCG | oligodT_307 CTGAGAGAATA | oligodT_319 AACCGGCCTT | oligodT_331 TTGTGAAGGC | oligodT_343 GCGAGATGGA | oligodT_355 AACAGTGACT | oligodT_367 AGAGGCACGA | oligodT_379 TGATTCCACA | 7 |
| 8 | oligodT_296 GATCAACAAG | oligodT_308 TGACAAGAGG | oligodT_320 CATTATCGCT | oligodT_332 GTCGACTCCT | oligodT_344 GTCATCGCGT | oligodT_356 CCACTTGGAT | oligodT_368 TTGGCTCCAA | oligodT_380 TAAGGGCATC | 8 |
| 9 | oligodT_297 ACCTTCGGAC | oligodT_309 TGTTTCGCTCT | oligodT_321 CATTGAGCTA | oligodT_333 CTGAGCTTGT | oligodT_345 AGTGGTAGCA | oligodT_357 TCAGACTGGT | oligodT_369 TTAAGTCCGT | oligodT_381 AATCCAGGAT | 9 |
| 10 | oligodT_298 TATCATGTGC | oligodT_310 AGGTATTCCGG | oligodT_322 CTTCGGAATC | oligodT_334 CCTCTGTCAT | oligodT_346 TATCCAGTTC | oligodT_358 GCATAGCCAA | oligodT_370 TCGGAGACCT | oligodT_382 AAGGTACCCT | 10 |
| 11 | oligodT_299 CGTCTAGTAA | oligodT_311 ACGGCCTACA | oligodT_323 AACCACACTA | oligodT_335 GCAGATGTAA | oligodT_347 AGTCGACGTA | oligodT_359 ACGTGCATCG | oligodT_371 TCCAGCGAAG | oligodT_383 ACGACAATGA | 11 |
| 12 | oligodT_300 TCTTAGTGTTG | oligodT_312 CACGTAAGAG | oligodT_324 GTACAGCGGA | oligodT_336 CTCATCATCT | oligodT_348 TTCGTTTACA | oligodT_360 GCCATGCAGA | oligodT_372 GTACTAAGAG | oligodT_384 TGCCTTCTGC | 12 |
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | |

Figure 9: Schematic representation of the distribution of barcoded oligodT primers in Array 4 of the cellenCHIP 384.

IV. Sequencing depth

The sequencing depth describes how many reads are obtained for each sample. In single cell RNA-sequencing this is normally scaled to the reads per single cell. Hence, the sequencing depth of a single library derived from the “cellenCHIP 384 – 3'RNA-seq Kit” requires 384 times the desired sequencing depth of each single cell. The sequencing depth cannot be adjusted differently for the different cells, hence each cell will roughly get the same amount of sequencing reads. Please note however, due to the early barcoding strategy of the “cellenCHIP 384 – 3'RNA-seq Kit” (i.e. pooling of all cDNAs before amplification), cells which had more mRNA in the beginning will also obtain more sequencing reads as their cDNA is overrepresented in the cDNA pool. This is particularly important when using highly heterogenous cell populations, e.g. PBMCs. In this case it is advised to plan with a higher sequencing depth per single cell as desired, so that also smaller cells will obtain the desired sequencing depth.

For scRNA-seq libraries in general, the sequencing depth can range between a few thousand up to a million reads per cell. It has been shown, that at 1 million reads per cell the saturation is reached with no additional information obtained when sequenced deeper. It should be noted however, that this relationship of sequencing depth and obtained information (e.g. the number of genes detected) is not linear but follows a saturation curve. Therefore, while the doubling from 10.000 to 20.000 reads might lead to an almost doubling of information, a doubling of the sequencing depth from 300.000 to 600.000 reads per single cell might only yield a few percent increase of information ([Figure 10](#)).

There is not a simple answer to the question of how deep to sequence. However, there are several considerations that should be taken into account.

- First and most obvious, the higher the sequencing depth the more expensive it is.
- Second, while a higher sequencing depth will result in more usable information so does a higher number of data points, i.e. cells.
- Lastly, it highly depends on the effect size of the investigated differences. For example, the difference between 2 completely different cell lines can be seen with a low number of cells and a low sequencing depth. The smaller the expected changes the more statistical power is required, and hence a higher sequencing depth and/or a higher number of cells is needed.

For the “cellenCHIP 384 – 3'RNA-seq Kit”, we have seen very good clustering performance when using a sequencing depth of ~50.000 reads per cell, however, to obtain reliable information for differential gene expression analysis we recommend to sequences about 150.000 reads per cell.

Please note, a single library can also be sequenced several times and the data of the different sequencing runs can be merged. This can allow for an initial lower sequencing depth and a further higher sequencing depth when analysis shows the requirement.




Figure 10: Number of detected genes per single cell (“Number of Genes”) in correlation to their corresponding raw sequencing reads (“Number of reads”). Each dot represents a single cell. The size of the dots represents the number of raw sequencing reads per single cell.